# Introcan Safety® 2 IV Catheter

## Multi Access Blood Control

**en** IV catheter designed to minimize inadvertent needle sticks and blood exposure. Made of PUR, radiopaque, with or without fixation wings.

For single use only

Do not sterilize

Sterilized using Ethylene Oxide.

Non-pyrogenic

Consult Electronic Instructions for Use. [www.bbraunusa.com](http://www.bbraunusa.com)

Use-by Date

Reorder number

Lot number

Single Sterile Barrier System


Single Sterile Barrier System with Protective Packaging Outside

Contents sterile in unopened, undamaged package. Not made with natural rubber latex, DEHP or PVC. Pictorials are for reference only. Pat.: [www.bbraunusa.com/patents](http://www.bbraunusa.com/patents)






**RADIOPAQUE**

Rx only

P-9223-1 REV. 12/22


15328481

**en**

Fig. 1: Before needle removal

- ① Flashback Chamber
- ② Protective Guard
- ③ Catheter Hub
- ④ Push-Off Plate
- ⑤ Needle Bevel
- ⑥ Wings
- ⑦ Safety Shield


Fig. 2: After needle removal

### en Description

The Introcan Safety® 2 IV Catheter is designed to:

- Minimize inadvertent needlestick injuries with a passive safety needle-shielding feature.
- Control the flow of blood which aids in the prevention of blood exposure.

The catheter may be used for any patient population with consideration given to adequacy of vascular anatomy and appropriateness for the solution being infused and duration of therapy.

#### Materials used

ABS, Chrome-Nickel Steel, PP, POM, PUR, Elastomer

#### Indication

- The Introcan Safety® 2 IV Catheter is inserted into a patient's vascular system for short term use to sample blood, monitor blood pressure, or administer fluids and blood intravascularly. The catheters may be used with power injectors at a maximum pressure of 325 psi with a luer lock connection only.

#### Contraindications

- Introcan Safety® 2 IV Catheter should not be used in patients with known hypersensitivity to any of the materials employed.

#### MRI Safety Information

#### MR Conditional

Non-clinical testing demonstrated that the Introcan Safety® 2 IV Catheter is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 T or 3.0 T
- Maximum spatial field gradient of 2,000 gauss/cm (20 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First level Controlled Operating Mode)

Under the scan conditions defined above, the Introcan Safety® 2 IV Catheter is expected to produce a maximum temperature rise of 2°C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Introcan Safety® 2 IV Catheter extends approximately 8 mm from this implant when imaged using a gradient echo pulse sequence and a 3.0 Tesla MR System.

#### Guidelines

- This device features a multi access blood control septum. After needle withdrawal, blood from catheter hub is restricted until connection of a luer access device which opens the fluid path again.
- Observe Standard Precautions on ALL patients. Aseptic technique, proper skin preparation, and continued protection of the insertion site are essential.
- Use of only ISO Luer-slip and Luer-lock connectors ensures proper functioning.
- This device is intended for single use only and is provided sterile and non-pyrogenic.
- This device is sterile unless packaging has been opened or damaged.

#### Risks & Warnings

- Re-use of single-use devices poses a potential risk to patient and user. It may lead to loss of function of the device, or contamination, which could cause severe injury or death.
- Always prime the catheter prior to arterial puncture.
- Verify adequate collateral circulation prior to arterial puncture.
- Non-ischemic complications of arterial catheterization, or even repeated attempts at catheterization, include bleeding, pseudoaneurysms, arteriovenous fistulae, nerve palsy, infection and injury to the tendon sheaths and adjacent structures during insertion.
- Arterial catheterization risks include injury to adjacent structures during insertion, such as nerve and tendon injury, infection, ischemic injury, such as vascular spasm and thromboembolism.
- Always clearly mark arterial lines to avoid accidental intra-arterial injection.
- Care must be taken to avoid needlesticks.
- Do not attempt to override or defeat the safety shield.
- In the unlikely event that the safety mechanism does not engage, keep the needle point away from the body and fingers at all times and immediately dispose into an approved sharps container.
- In the case of an unsuccessful IV start, remove the needle first to activate the safety mechanism, then remove catheter from patient and discard both.
- Report needlestick injuries immediately and follow established institutional protocols.
- Do not bend the catheter or needle during insertion, advancement, or removal of the needle.
- Extreme care should be taken not to cut the catheter and possibly cause embolization.
- Do not use scissors or sharp instruments at or near the insertion site.
- Never reinsert the needle inside the catheter once the needle has been partially or completely withdrawn as it may pierce and/or sever the catheter.
- Luer-slip connections should not be left unattended due to potential of disconnection.
- If blood is present inside the catheter, flush the device according to institutional protocol.
- Not intended to facilitate or aid in the placement of vascular access devices such as guidewires, indwelling central venous catheters (CVC), peripherally inserted central catheters (PICC) and midline catheters into the vascular system.

#### Pressure Injection Information/Warnings and Cautions

- Prior to use power injectors, ensure that a secure connection exists between the catheter and power injector.
- Patency of the catheter must be assured immediately before power injecting.
- Measures should be taken to avoid kinking or obstruction of the catheter system during power injection to avoid product failure.

- Pressure limit settings on high pressure injector equipment may not prevent over pressurizing an occluded or partially occluded catheter, which may cause catheter failure.
- Follow the contrast media manufactures instructions for use, contraindications, warnings and precautions.
- Due to variations in add-on devices, tubing, injectate, temperature and pressure limit settings, the listed injection flow rates may not be achievable.

#### Power Injection Flow Rates

Gauge / Length	Flow Rate (mL/s) Contrast Media Viscosity at 20°C: 2.3 mPa*s	Flow Rate (mL/s) Contrast Media Viscosity at 20°C: 27.5 mPa*s
24 Ga. x 0.55 in.	5.0	2.5
24 Ga. x ¾ in.	5.0	2.5
22 Ga. x 1 in.	7.0	3.5
20 Ga. x 1 in.	9.5	4.0
20 Ga. x 1 ¼ in.	9.5	4.0
20 Ga. x 2 in.	9.5	4.0
18 Ga. x 1 ¼ in.	12.5	5.0
18 Ga. x 1 ¾ in.	12.5	5.0

#### Duration of Use

- Change device according to CDC Guidelines and/or institutional protocol.
- The insertion site should be checked frequently and at regular intervals.
- Introcan Safety® 2 IV Catheter should be removed in the event of local or systemic signs of infection.

#### Instructions for Use: Use aseptic technique. Do not rotate catheter hub prior to insertion.

1. Gather the necessary supplies.
2. Select and prepare site per institutional protocol.
3. Apply a tourniquet.
4. Holding flashback chamber, remove protective guard in a straight outward motion. Inspect device.
5. Hold the device by the flashback chamber with thumb and fingers on opposite sides. Verify push-off plate and needle bevel are in the "up" position.
6. Anchor the vessel with gentle skin traction. Adjust angle of insertion and puncture a suitable vessel. If vessel puncture was successful, blood will immediately be visible inside the transparent flashback chamber (see figure A).
7. Lower and advance the entire catheter and needle unit slightly (approximately 1/8 in. or 3 mm) to ensure catheter tip is in the vessel (see figure B).
8. Using the push-off plate, advance the catheter off the needle and observe for blood return between the needle and catheter (see figure C). After confirmation, continue advancing catheter off the needle and into the vessel.
9. Remove tourniquet.
10. Stabilize the catheter hub to prevent catheter dislodgement during needle removal. Withdraw the needle straight out with a controlled and continuous motion (minimize rotation or bending of the needle). The metal safety shield will automatically attach to the needle tip as the needle exits the catheter hub (see figure D). Blood flow from catheter hub is restricted after needle is removed and until connection of a Luer access device is made.
11. Dispose of the shielded needle immediately into an approved sharps container.
12. Immediately connect the infusion line or accessory device and cover puncture site with a sterile and transparent dressing (see figure E) per institutional protocol.
13. After removal, dispose of the catheter per institutional protocol.

Distributed by:  
**B. Braun Medical Inc.**  
Bethlehem, PA 18018-3524 USA  
1-800-227-2862

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**B. Braun Melsungen AG**